

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Monkey CYP2R1 Polyclonal Antibody

多克隆抗体

产品货号: PA03400

第三版

**描述:** This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This enzyme is a microsomal vitamin D hydroxylase that converts vitamin D into the active ligand for the vitamin D receptor. A mutation in this gene has been associated with selective 25-hydroxyvitamin D deficiency. CYP2R1 (Cytochrome P450 Family 2 Subfamily R Member 1) is a Protein Coding gene. Diseases associated with CYP2R1 include Rickets Due To Defect In Vitamin D 25-Hydroxylation and Hypocalcemic Vitamin D-Dependent Rickets. Among its related pathways are Cytochrome P450 - arranged by substrate type and cholesterol biosynthesis I. GO annotations related to this gene include iron ion binding and oxidoreductase activity, acting on paired donors, with incorporation or reduction of molecular oxygen. An important paralog of this gene is CYP2J2.

**抗原:** Synthesized peptide derived from the Internal region of human CYP2R1

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) , 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:**