

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human /Mouse /Rat Gabra3 Polyclonal Antibody

多克隆抗体

产品货号: PA03417

第三版

描述:GABA is the major inhibitory neurotransmitter in the mammalian brain where it acts at GABA-A receptors, which are ligand-gated chloride channels. Chloride conductance of these channels can be modulated by agents such as benzodiazepines that bind to the GABA-A receptor. At least 16 distinct subunits of GABA-A receptors have been identified. GABRA3 (Gamma-Aminobutyric Acid Type A Receptor Alpha3 Subunit) is a Protein Coding gene. Diseases associated with GABRA3 include Thyrotoxic Periodic Paralysis and Mucopolysaccharidosis II. Among its related pathways are Ligand-gated ion channel transport and Transport of glucose and other sugars, bile salts and organic acids, metal ions and amine compounds. GO annotations related to this gene include chloride channel activity and GABA-A receptor activity. An important paralog of this gene is GABRA1.

抗原:Synthetic peptide of Rat Gabra3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: