
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CK-10 Polyclonal Antibody

多克隆抗体

产品货号: PA03458

第三版

描述: Cytokeratin 10 (CK10) has a molecular weight of 56.5KDa and is mainly marked on the basal and granular cell layers of the epithelium. Meanwhile, the expression of CK10 is proportional to the degree of cell differentiation, and the highly differentiated ones are often positive, which is mainly used in the diagnosis of squamous cell carcinoma.

抗原: Synthesized peptide derived from the Internal region of human Cytokeratin 10

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Seen in all suprabasal cell layers including stratum corneum.