

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CK-13 Polyclonal Antibody

多克隆抗体

产品货号: PA03459

第三版

描述:Cytokeratin 13 (CK13) has a molecular weight of 53 KDa, it is mainly marked on non-keratinized squamous epithelium and can be used for diagnosis of squamous cell carcinoma.

抗原:Synthesized peptide derived from the Internal region of human Cytokeratin 13.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in some epidermal sweat gland ducts (at protein level) and in exocervix, esophagus and placenta.