

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human /Mouse TFDP1 Polyclonal Antibody

多克隆抗体

产品货号: PA03475

第三版

描述: This gene encodes a member of a family of transcription factors that heterodimerize with E2F proteins to enhance their DNA-binding activity and promote transcription from E2F target genes. The encoded protein functions as part of this complex to control the transcriptional activity of numerous genes involved in cell cycle progression from G1 to S phase. Alternative splicing results in multiple transcript variants. Pseudogenes of this gene are found on chromosomes 1, 15, and X. TFDP1 (Transcription Factor Dp-1) is a Protein Coding gene. Diseases associated with TFDP1 include Retinoblastoma. Among its related pathways are Cellular Senescence and Transcriptional activity of SMAD2/SMAD3-SMAD4 heterotrimer. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and transcription coactivator activity. An important paralog of this gene is TFDP3.

抗原: Synthetic peptide of human TFDP1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Highest levels in muscle. Also expressed in brain, placenta, liver and kidney. Lower levels in lung and pancreas. Not detected in heart.