
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Monkey DBF4 Polyclonal Antibody

多克隆抗体

产品货号: PA03500

第三版

描述:DBF4 (DBF4 Zinc Finger) is a Protein Coding gene. Among its related pathways are E2F mediated regulation of DNA replication and Cell cycle. GO annotations related to this gene include nucleic acid binding and enzyme activator activity. An important paralog of this gene is DBF4B.

抗原:Synthesized peptide derived from the N-terminal region of human Dbf4

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Highly expressed in testis and thymus. Expressed also in most cancer cells lines.