

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse DCLK2 Polyclonal Antibody

多克隆抗体

产品货号: PA03503

第三版

描述:Protein kinase with a significantly reduced Ca²⁺/CAM affinity and dependence compared to other members of the CaMK family. May play a role in the down-regulation of CRE-dependent gene activation probably by phosphorylation of the CREB coactivator CRT2/TORC2 and the resulting retention of TORC2 in the cytoplasm.

抗原:Synthesized peptide derived from the N-terminal region of human DCAMKL2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in the brain, heart and eyes.