

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human OR51E1 Polyclonal Antibody

多克隆抗体

产品货号: PA03520

第三版

描述:Odorant receptor.

抗原:Synthesized peptide derived from the C-terminal region of human D-GPCR

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Highly expressed in prostate. Very low levels may be detected in some other tissues, such as placenta, skeletal muscle, heart, ovary and testis. Up-regulated in prostate cancers.