
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ZDHHC15 Polyclonal Antibody

多克隆抗体

产品货号: PA03521

第三版

描述:The protein encoded by this gene belongs to the DHHC palmitoyltransferase family. Mutations in this gene are associated with mental retardatio X-linked type 91 (MRX91). Alternatively spliced transcript variants encoding different isoforms have been found for this gene. ZDHHC15 (Zinc Finger DHHC-Type Containing 15) is a Protein Coding gene. GO annotations related to this gene include receptor activity and palmitoyltransferase activity. An important paralog of this gene is ZDHHC2.

抗原:Synthesized peptide derived from the C-terminal region of human DHHC-15

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in placenta, liver, lung, kidney, heart and brain.