

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human DHODH Polyclonal Antibody

多克隆抗体

产品货号: PA03524

第三版

描述:The protein encoded by this gene catalyzes the fourth enzymatic step, the ubiquinone-mediated oxidation of dihydroorotate to orotate, in de novo pyrimidine biosynthesis. This protein is a mitochondrial protein located on the outer surface of the inner mitochondrial membrane.

抗原:Synthesized peptide derived from the Internal region of human DHODH

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: