

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat DLX5 Polyclonal Antibody

多克隆抗体

产品货号: PA03533

第三版

描述: This gene encodes a member of a homeobox transcription factor gene family similar to the Drosophila distal-less gene. The encoded protein may play a role in bone development and fracture healing. Mutation in this gene, which is located in a tail-to-tail configuration with another member of the family on the long arm of chromosome 7, may be associated with split-hand/split-foot malformation. DLX5 (Distal-Less Homeobox 5) is a Protein Coding gene. Diseases associated with DLX5 include Split-Hand/Foot Malformation 1 With Sensorineural Hearing Loss and Split Hand Foot Malformation 1. Among its related pathways are Signaling pathways regulating pluripotency of stem cells and Neural Crest Differentiation. GO annotations related to this gene include RNA polymerase II core promoter proximal region sequence-specific DNA binding and transcriptional activator activity, RNA polymerase II core promoter proximal region sequence-specific binding. An important paralog of this gene is DLX3.

抗原: Synthesized peptide derived from the Internal region of human Dlx-5

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: