

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human /Mouse /Rat CUEDC2 Polyclonal Antibody

多克隆抗体

产品货号: PA03543

第三版

**描述:**CUEDC2 is involved in ubiquitin- and proteasome-mediated degradation of progesterone receptor and estrogen receptor (ER)-alpha. CUEDC2 reduced PRB protein content and promoted progesterone-induced PRB degradation via the ubiquitin-proteasome pathway.

**抗原:**Synthetic peptide of human CUEDC2

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Significantly up-regulated in breast tumor tissues compared with matched adjacent normal tissues (at protein level). Levels inversely correlate with ESR1 in breast cancers and are lower in low-grade tumors compared to high-grade tumors.