
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse DOK4 Polyclonal Antibody

多克隆抗体

产品货号: PA03572

第三版

描述:DOK4 (docking protein 4) belongs to the family of docking proteins which interact with receptor tyrosine kinases and mediate particular biological responses. DOK4 functions in RET mediated neurite outgrowth and plays a positive role in activation of the MAP kinase pathway.

抗原:Synthesized peptide derived from the N-terminal region of human Dok-4

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Widely expressed. High expression in skeletal muscle, heart, kidney and liver. Weaker expression in spleen, lung and small intestine, brain, heart and. Expressed in both resting and activated peripheral blood T-cells.