

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ZFYVE16 Polyclonal Antibody

多克隆抗体

产品货号: PA03590

第三版

描述: This gene encodes an endosomal protein that belongs to the FYVE zinc finger family of proteins. The encoded protein is thought to regulate membrane trafficking in the endosome. This protein functions as a scaffold protein in the transforming growth factor-beta signaling pathway and is involved in positive and negative feedback regulation of the bone morphogenetic protein signaling pathway. Alternate splicing results in multiple transcript variants.

抗原: Synthetic peptide of human ZFYVE16

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Widely expressed. Highly expressed in kidney, placenta and lung. Expressed at intermediate level in heart, brain, skeletal muscle, spleen and liver. Weakly expressed in colon, thymus and peripheral blood lymphocytes.