

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human /Mouse ZHX3 Polyclonal Antibody

多克隆抗体

产品货号: PA03600

第三版

描述: This gene encodes a member of the zinc fingers and homeoboxes (ZHX) gene family. The encoded protein contains two C2H2-type zinc fingers and five homeodomains and forms a dimer with itself or with zinc fingers and homeoboxes family member 1. In the nucleus, the dimerized protein interacts with the A subunit of the ubiquitous transcription factor nuclear factor- κ B and may function as a transcriptional repressor.

抗原: Synthetic peptide of human ZHX3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20°C 保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4)，1% BSA

应用: WB 1 ~ 5 $\mu\text{g/ml}$.

特异性: Widely expressed. High expression in kidney.