



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human DTYMK Polyclonal Antibody

多克隆抗体

产品货号: PA03612

第三版

**描述:**DTYMK (Deoxythymidylate Kinase) is a Protein Coding gene. Among its related pathways are superpathway of pyrimidine deoxyribonucleotides de novo biosynthesis and Purine metabolism (KEGG). GO annotations related to this gene include kinase activity and thymidylate kinase activity.

**抗原:**Synthesized peptide derived from the N-terminal region of human dTMP Kinase

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**