

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human EAAT1 Polyclonal Antibody

多克隆抗体

产品货号: PA03632

第三版

**描述:** This gene encodes a member of a high affinity glutamate transporter family. This gene functions in the termination of excitatory neurotransmission in central nervous system. Mutations are associated with episodic ataxia, Type 6. Alternative splicing results in multiple transcript variants.

**抗原:** Synthesized peptide derived from the C-terminal region of human EAAT1.

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Highly expressed in cerebellum, but also found in frontal cortex, hippocampus and basal ganglia.