

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat EAR2 Polyclonal Antibody

多克隆抗体

产品货号: PA03634

第三版

描述:NR2F6 (Nuclear Receptor Subfamily 2 Group F Member 6) is a Protein Coding gene. Among its related pathways are Circadian rhythm related genes and Gene Expression. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and RNA polymerase II distal enhancer sequence-specific DNA binding. An important paralog of this gene is NR2F1.

抗原:Synthesized peptide derived from the N-terminal region of human EAR2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in heart, placenta, liver, skeletal muscle, kidney and pancreas.