

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human EGFL4 Polyclonal Antibody

多克隆抗体

产品货号: PA03669

第三版

描述:The protein encoded by this gene is a single-pass type I membrane protein of unknown function that contains several EGF-like domains, Kelch repeats, and PSI domains. Defects in this gene are a cause of Carpenter syndrome 2. Two transcript variants encoding different isoforms have been found for this gene. MEGF8 (Multiple EGF Like Domains 8) is a Protein Coding gene. Diseases associated with MEGF8 include Carpenter Syndrome 2 and Carpenter Syndrome. GO annotations related to this gene include calcium ion binding and receptor activity. An important paralog of this gene is ATRNL1.

抗原:Synthesized peptide derived from the Internal region of human EGFL4

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: