
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Mouse/Rat PAFAH1B2 Polyclonal Antibody

多克隆抗体

产品货号: PA03721

第三版

描述:Platelet-activating factor acetylhydrolase (PAFAH) inactivates platelet-activating factor (PAF) into acetate and LYSO-PAF. This gene encodes the beta subunit of PAFAH, the other subunits are alpha and gamma. Multiple alternatively spliced transcript variants have been described for this gene.

抗原:Recombinant human Platelet activating factor acetylhydrolase 1b catalytic subunit 2 protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitous.