

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Mouse THRAP3 Polyclonal Antibody

多克隆抗体

产品货号: PA03740

第三版

描述:THRAP3 (Thyroid Hormone Receptor Associated Protein 3) is a Protein Coding gene. Among its related pathways are PEDF Induced Signaling and Circadian rythm related genes. GO annotations related to this gene include poly(A) RNA binding and transcription coactivator activity. An important paralog of this gene is BCLAF1.

抗原:Recombinant human Thyroid hormone receptor-associated protein 3 protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitous.