

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human EXOG Polyclonal Antibody

多克隆抗体

产品货号: PA03748

第三版

描述:This gene encodes an endo/exonuclease with 5'-3' exonuclease activity. The encoded enzyme catalyzes the hydrolysis of ester linkages at the 5' end of a nucleic acid chain. This enzyme is localized to the mitochondria and may play a role in programmed cell death. Alternatively spliced transcript variants have been described. A pseudogene exists on chromosome 18.EXOG (Exo/Endonuclease G) is a Protein Coding gene. GO annotations related to this gene include nucleic acid binding and endonuclease activity. An important paralog of this gene is ENDOG.

抗原:Synthesized peptide derived from the Internal region of human Endo G-L1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitous.