

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Mouse/Rat EIF5A2 Polyclonal Antibody

多克隆抗体

产品货号: PA03762

第三版

描述:EIF5A2 (Eukaryotic Translation Initiation Factor 5A2) is a Protein Coding gene. Among its related pathways are Gamma carboxylation, hypusine formation and arylsulfatase activation and Metabolism of proteins. GO annotations related to this gene include RNA binding and translation elongation factor activity. An important paralog of this gene is EIF5A.

抗原:Recombinant human Eukaryotic translation initiation factor 5A-2 protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in ovarian and colorectal cancer cell lines (at protein level). Highly expressed in testis. Overexpressed in some cancer cells.