
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CAMK4 Polyclonal Antibody

多克隆抗体

产品货号: PA03766

第三版

描述:The product of this gene belongs to the serine/threonine protein kinase family, and to the Ca(2+)/calmodulin-dependent protein kinase subfamily. This enzyme is a multifunctional serine/threonine protein kinase with limited tissue distribution, that has been implicated in transcriptional regulation in lymphocytes, neurons and male germ cells.

抗原:Recombinant human Calcium/calmodulin-dependent protein kinase type IV protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in epithelial ovarian cancer tissue.