

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse EPS8L1 Polyclonal Antibody

多克隆抗体

产品货号: PA03800

第三版

描述: This gene encodes a protein that is related to epidermal growth factor receptor pathway substrate 8 (EPS8), a substrate for the epidermal growth factor receptor. The function of this protein is unknown. At least two alternatively spliced transcript variants encoding different isoforms have been found for this gene. EPS8L1 (EPS8 Like 1) is a Protein Coding gene. GO annotations related to this gene include actin binding and Rho guanyl-nucleotide exchange factor activity. An important paralog of this gene is EPS8.

抗原: Synthesized peptide derived from the Internal region of human Eps8L1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4)，1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Detected in placenta.