

---

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse EPS8L2 Polyclonal Antibody

多克隆抗体

产品货号: PA03805

---

第三版

**描述:**This gene encodes a member of the EPS8 gene family. The encoded protein, like other members of the family, is thought to link growth factor stimulation to actin organization, generating functional redundancy in the pathways that regulate actin cytoskeletal remodeling. EPS8L2 (EPS8 Like 2) is a Protein Coding gene. GO annotations related to this gene include actin binding and Rac guanyl-nucleotide exchange factor activity. An important paralog of this gene is EPS8.

**抗原:**Synthesized peptide derived from the N-terminal region of human Eps8L2

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4)，1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**Detected in fibroblasts and placenta.