

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat EPN1 Polyclonal Antibody

多克隆抗体

产品货号: PA03806

第三版

**描述:** This gene encodes a member of the epsin protein family. The encoded protein binds clathrin and is involved in the endocytosis of clathrin-coated vesicles. Loss of function of this gene is associated with reduced tumor growth and progression in certain cancer types. EPN1 (Epsin 1) is a Protein Coding gene. Diseases associated with EPN1 include Erythrocytosis, Somatic. Among its related pathways are Sertoli-Sertoli Cell Junction Dynamics and Notch signaling pathway (KEGG). GO annotations related to this gene include transcription factor binding and phospholipid binding. An important paralog of this gene is EPN2.

**抗原:** Synthesized peptide derived from the Internal region of human Epsin 1

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**