

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat ERK 3 Polyclonal Antibody

多克隆抗体

产品货号: PA03824

---

第三版

**描述:** Phosphorylates microtubule-associated protein 2 (MAP2). May promote entry in the cell cycle.

**抗原:** Synthesized peptide derived from the Internal region of human ERK 3

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Highest expression in the skeletal muscle, followed by the brain. Also found in heart, placenta, lung, liver, pancreas, kidney and skin fibroblasts.