

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat ERK 3 Polyclonal Antibody

多克隆抗体

产品货号: PA03840

第三版

描述: Phosphorylates microtubule-associated protein 2 (MAP2). May promote entry in the cell cycle.

抗原: Synthesized peptide derived from human ERK 3 around the non-phosphorylation site of Ser189.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Highest expression in the skeletal muscle, followed by the brain. Also found in heart, placenta, lung, liver, pancreas, kidney and skin fibroblasts.