

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat TRIM72 Polyclonal Antibody

多克隆抗体

产品货号: PA03880

第三版

描述:TRIM72 (Tripartite Motif Containing 72) is a Protein Coding gene. Among its related pathways are Cardiac conduction and Smooth Muscle Contraction. GO annotations related to this gene include phosphatidylserine binding. An important paralog of this gene is TRIM50.

抗原:Recombinant Protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: