

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat MAP1LC3B Polyclonal Antibody

多克隆抗体

产品货号: PA03893

第三版

描述:The product of this gene is a subunit of neuronal microtubule-associated MAP1A and MAP1B proteins, which are involved in microtubule assembly and important for neurogenesis. Studies on the rat homolog implicate a role for this gene in autophagy, a process that involves the bulk degradation of cytoplasmic component. MAP1LC3B (Microtubule Associated Protein 1 Light Chain 3 Beta) is a Protein Coding gene. Diseases associated with MAP1LC3B include Sporadic Pheochromocytoma and Lacrimal Gland Adenocarcinoma. Among its related pathways are Vesicle-mediated transport and Senescence and Autophagy in Cancer. GO annotations related to this gene include microtubule binding. An important paralog of this gene is MAP1LC3B2.

抗原:Recombinant Protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Most abundant in heart, brain, skeletal muscle and testis. Little expression observed in liver.