

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Mouse/Rat CACNB3 Polyclonal Antibody

多克隆抗体

产品货号: PA03899

第三版

描述: This gene encodes a regulatory beta subunit of the voltage-dependent calcium channel. Beta subunits are composed of five domains, which contribute to the regulation of surface expression and gating of calcium channels and may also play a role in the regulation of transcription factors and calcium transport. CACNB3 (Calcium Voltage-Gated Channel Auxiliary Subunit Beta 3) is a Protein Coding gene. Among its related pathways are CCR5 Pathway in Macrophages and PEDF Induced Signaling. GO annotations related to this gene include protein kinase binding and high voltage-gated calcium channel activity. An important paralog of this gene is CACNB4.

抗原: Synthetic Peptide

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed mostly in brain, smooth muscle and ovary.