

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human FASN Polyclonal Antibody

多克隆抗体

产品货号: PA03905

第三版

描述:The enzyme encoded by this gene is a multifunctional protein. Its main function is to catalyze the synthesis of palmitate from acetyl-CoA and malonyl-CoA, in the presence of NADPH, into long-chain saturated fatty acids. In some cancer cell lines, this protein has been found to be fused with estrogen receptor-alpha (ER-alpha), in which the N-terminus of FAS is fused in-frame with the C-terminus of ER-alpha.

抗原:Synthesized peptide derived from the Internal region of human Fatty Acid Synthase.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitous. Prominent expression in brain, lung, and liver.