

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human FGL2 Polyclonal Antibody

多克隆抗体

产品货号: PA03925

第三版

描述:The protein encoded by this gene is a secreted protein that is similar to the beta- and gamma-chains of fibrinogen. The carboxyl-terminus of the encoded protein consists of the fibrinogen-related domains (FRED). The encoded protein forms a tetrameric complex which is stabilized by interchain disulfide bonds. This protein may play a role in physiologic functions at mucosal sites.FGL2 (Fibrinogen Like 2) is a Protein Coding gene. Diseases associated with FGL2 include Viral Hepatitis and Hepatitis B. Among its related pathways are Innate Immune System. An important paralog of this gene is ANGPT2.

抗原:Synthesized peptide derived from the N-terminal region of human Fgl2.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Constitutively expressed in cytotoxic T-cells.