

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Fibulin 2 Polyclonal Antibody

多克隆抗体

产品货号: PA03927

第三版

描述: This gene encodes an extracellular matrix protein, which belongs to the fibulin family. This protein binds various extracellular ligands and calcium. It may play a role during organ development, in particular, during the differentiation of heart, skeletal and neuronal structures. Alternatively spliced transcript variants encoding different isoforms have been identified.

抗原: Synthesized peptide derived from the Internal region of human Fibulin-2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Component of both basement membranes and other connective tissues. Expressed in heart, placenta and ovary.