

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Fibulin 3 Polyclonal Antibody

多克隆抗体

产品货号: PA03928

第三版

描述: This gene encodes a member of the fibulin family of extracellular matrix glycoproteins. Like all members of this family, the encoded protein contains tandemly repeated epidermal growth factor-like repeats followed by a C-terminus fibulin-type domain. This gene is upregulated in malignant gliomas and may play a role in the aggressive nature of these tumors. Mutations in this gene are associated with Doyme honeycomb retinal dystrophy. Alternatively spliced transcript variants that encode the same protein have been described. EFEMP1 (EGF Containing Fibulin Like Extracellular Matrix Protein 1) is a Protein Coding gene. Diseases associated with EFEMP1 include Doyme Honeycomb Degeneration Of Retina and Familial Drusen. Among its related pathways are Phospholipase-C Pathway and ERK Signaling. GO annotations related to this gene include calcium ion binding and epidermal growth factor receptor binding. An important paralog of this gene is EFEMP2.

抗原: Synthesized peptide derived from the Internal region of human Fibulin-3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: In the eye, associated with photoreceptor outer and inner segment regions, the nerve fiber layer, outer nuclear layer and inner and outer plexiform layers of the retina.