

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse FARP2 Polyclonal Antibody

多克隆抗体

产品货号: PA03931

第三版

**描述:**FARP2 (FERM, ARH/RhoGEF And Pleckstrin Domain Protein 2) is a Protein Coding gene. Diseases associated with FARP2 include Pseudohypoparathyroidism Ia and Aarskog-Scott Syndrome. Among its related pathways are Development Endothelin-1/EDNRA signaling and Focal Adhesion. GO annotations related to this gene include Rho guanyl-nucleotide exchange factor activity and Rac guanyl-nucleotide exchange factor activity. An important paralog of this gene is FARP1.

**抗原:**Synthesized peptide derived from the Internal region of human FIR

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**