

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse TNS2 Polyclonal Antibody

多克隆抗体

产品货号: PA03963

第三版

描述:The protein encoded by this gene belongs to the tensin family. Tensin is a focal adhesion molecule that binds to actin filaments and participates in signaling pathways. This protein plays a role in regulating cell migration. Alternative splicing occurs at this locus and three transcript variants encoding three distinct isoforms have been identified. TNS2 (Tensin 2) is a Protein Coding gene. Among its related pathways are Adhesion and PI3K / Akt Signaling. GO annotations related to this gene include kinase binding and phosphoprotein phosphatase activity. An important paralog of this gene is TNS1.

抗原:Synthesized peptide derived from human Tensin-2 around the non-phosphorylation site of Tyr483.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Detected in heart, kidney, brain, thymus, spleen, liver, placenta, lung, skeletal muscle and small intestine.