

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse FOXB1/2 Polyclonal Antibody

多克隆抗体

产品货号: PA03968

第三版

**描述:**FOXB1 is a winged helix/forkhead transcription factor. FOXB1 is specifically expressed in the developing central nervous system (CNS). Early embryonic FOXB1 expression is restricted to the mammillary body region of the caudal hypothalamus, midbrain, hindbrain and spinal cord. FOXB1 may play a role in postnatal growth, lactation and CNS development.

**抗原:**Synthesized peptide derived from the Internal region of human FoxB1/2

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**