

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse MOB3A Polyclonal Antibody

多克隆抗体

产品货号: PA03973

第三版

描述: MOB3A (MOB Kinase Activator 3A) is a Protein Coding gene. An important paralog of this gene is MOB3B. MOB3B (MOB Kinase Activator 3B) is a Protein Coding gene. An important paralog of this gene is MOB3A. The protein encoded by this gene shares similarity with the yeast Mob1 protein. Yeast Mob1 binds Mps1p, a protein kinase essential for spindle pole body duplication and mitotic checkpoint regulation. This gene is located on the opposite strand as the interferon kappa precursor (IFNK) gene.

抗原: Synthesized peptide derived from the N-terminal region of human Mob3A/B

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4)，1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: