
本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat HTR1D Polyclonal Antibody

多克隆抗体

产品货号: PA03974

第三版

描述:HTR1D (5-Hydroxytryptamine Receptor 1D) is a Protein Coding gene. Diseases associated with HTR1D include Migraine With Or Without Aura 1 and Obsessive-Compulsive Disorder. Among its related pathways are Serotonergic synapse and Peptide ligand-binding receptors. GO annotations related to this gene include G-protein coupled receptor activity and serotonin binding. An important paralog of this gene is HTR1B.

抗原:Synthesized peptide derived from the Internal region of human SR-1D

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Detected in brain neocortex and caudate nucleus (at protein level).