
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL1F10 Polyclonal Antibody

多克隆抗体

产品货号: PA03993

第三版

描述:IL1F10 is a member of the interleukin 1 cytokine family. This cytokine is thought to participate in a network of interleukin 1 family members to regulate adapted and innate immune responses.

抗原:Synthesized peptide derived from the C-terminal region of human IL-1F10.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in fetal skin, spleen and tonsil. Expressed mostly in the basal epithelia of skin and in proliferating B-cells of the tonsil.