

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CMTM8 Polyclonal Antibody

多克隆抗体

产品货号: PA03999

第三版

描述:This gene belongs to the chemokine-like factor gene superfamily, a novel family that is similar to the chemokine and the transmembrane 4 superfamilies. This gene is one of several chemokine-like factor genes located in a cluster on chromosome 3. This gene is widely expressed in many tissues, but the exact function of the encoded protein is unknown.

抗原:Synthesized peptide derived from the Internal region of human CMTM8.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Highly expressed in liver and pancreas.