

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TP53INP2 Polyclonal Antibody

多克隆抗体

产品货号: PA04000

第三版

**描述:**The protein encoded by this gene promotes autophagy and is essential for proper autophagosome formation and processing. In addition, the encoded protein can enhance rDNA transcription by helping in the assembly of the POLR1/RNA polymerase I preinitiation complex. Finally, this protein serves as a transcriptional activator for some genes. TP53INP2 (Tumor Protein P53 Inducible Nuclear Protein 2) is a Protein Coding gene. Diseases associated with TP53INP2 include Obesity. Among its related pathways are Autophagy - animal. An important paralog of this gene is TP53INP1.

**抗原:**Synthesized peptide derived from the C-terminal region of human TP53INP2.

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μ g/ml.

**特异性:**