

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat TMPO Polyclonal Antibody

多克隆抗体

产品货号: PA04016

第三版

描述:TMPO (Thymopoietin) is a Protein Coding gene. Diseases associated with TMPO include Cardiomyopathy, Dilated, 1T and Tmpo-Related Dilated Cardiomyopathy. Among its related pathways are Nuclear Envelope Reassembly and Apoptosis and Autophagy. GO annotations related to this gene include lamin binding. The protein encoded by this gene resides in the nucleus and may play a role in the assembly of the nuclear lamina, and thus help maintain the structural organization of the nuclear envelope. It may function as a receptor for the attachment of lamin filaments to the inner nuclear membrane. Mutations in this gene are associated with dilated cardiomyopathy. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

抗原:Synthetic peptide from human protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in many tissues