

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CSF2 Polyclonal Antibody

多克隆抗体

产品货号: PA04020

第三版

描述:CSF2 (Colony Stimulating Factor 2) is a Protein Coding gene. Diseases associated with CSF2 include Mucositis and Candida Glabrata. Among its related pathways are RET signaling and IL-15 Signaling Pathways and their Primary Biological Effects in Different Immune Cell Types. GO annotations related to this gene include cytokine activity and granulocyte macrophage colony-stimulating factor receptor binding. The protein encoded by this gene is a cytokine that controls the production, differentiation, and function of granulocytes and macrophages. The active form of the protein is found extracellularly as a homodimer. This gene has been localized to a cluster of related genes at chromosome region 5q31, which is known to be associated with interstitial deletions in the 5q- syndrome and acute myelogenous leukemia. Other genes in the cluster include those encoding interleukins 4, 5, and 13.

抗原:Synthetic peptide from human protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: