

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MZB1 Polyclonal Antibody

多克隆抗体

产品货号: PA04030

第三版

描述:ADCYAP1 (Adenylate Cyclase Activating Polypeptide 1) is a Protein Coding gene. Diseases associated with ADCYAP1 include Persian Gulf Syndrome and Holoprosencephaly-3. Among its related pathways are Presynaptic function of Kainate receptors and Sudden Infant Death Syndrome (SIDS) Susceptibility Pathways. GO annotations related to this gene include receptor binding and signal transducer activity, downstream of receptor. An important paralog of this gene is VIP. This gene encodes a secreted proprotein that is further processed into multiple mature peptides. These peptides stimulate adenylate cyclase and increase cyclic adenosine monophosphate (cAMP) levels, resulting in the transcriptional activation of target genes. The products of this gene are key mediators of neuroendocrine stress responses. Alternative splicing results in multiple transcript variants.

抗原:Synthetic peptide from human protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Widely expressed with highest levels in adult brain, small intestine and lymphoid tissues such as thymus and spleen. Expression is frequently lower in intestinal-type gastric cancer. In obese patients, more abundant in omental than in subcutaneous fat.