

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse GLI1 Polyclonal Antibody

多克隆抗体

产品货号: PA04034

第三版

描述:GLI1 (GLI Family Zinc Finger 1) is a Protein Coding gene. Diseases associated with GLI1 include Basal Cell Carcinoma and Necrobiosis Lipoidica. Among its related pathways are Wnt Signaling Pathways: beta-Catenin-independent Wnt/Ca²⁺ Signaling and Other Non-canonical Wnt Signaling Pathways and G-Beta Gamma Signaling. GO annotations related to this gene include chromatin binding and microtubule binding. An important paralog of this gene is GLI2. This gene encodes a member of the Kruppel family of zinc finger proteins. The encoded transcription factor is activated by the sonic hedgehog signal transduction cascade and regulates stem cell proliferation. The activity and nuclear localization of this protein is negatively regulated by p53 in an inhibitory loop. Multiple transcript variants encoding different isoforms have been found for this gene.

抗原:Synthesized peptide derived from human GLI-1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Testis, myometrium and fallopian tube. Also expressed in the brain with highest expression in the cerebellum, optic nerve and olfactory tract.