

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Acetyl-ACSS1 (Lys642) Polyclonal Antibody

多克隆抗体

产品货号: PA04035

第三版

描述: This gene encodes a mitochondrial acetyl-CoA synthetase enzyme. A similar protein in mice plays an important role in the tricarboxylic acid cycle by catalyzing the conversion of acetate to acetyl CoA. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. ACSS1 (Acyl-CoA Synthetase Short-Chain Family Member 1) is a Protein Coding gene. Among its related pathways are Statin Pathway and Sterol Regulatory Element-Binding Proteins (SREBP) signalling. GO annotations related to this gene include AMP binding and acetate-CoA ligase activity. An important paralog of this gene is ACSS2.

抗原: Synthesized Acetyl peptide derived from human ACSS1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: