

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ATG16L1 Polyclonal Antibody

多克隆抗体

产品货号: PA04060

第三版

描述:ATG16L1 (Autophagy Related 16 Like 1) is a Protein Coding gene. Diseases associated with ATG16L1 include Inflammatory Bowel Disease 10 and Inflammatory Bowel Disease. Among its related pathways are Autophagy Pathway and Senescence and Autophagy in Cancer. GO annotations related to this gene include identical protein binding. An important paralog of this gene is ATG16L2. The protein encoded by this gene is part of a large protein complex that is necessary for autophagy, the major process by which intracellular components are targeted to lysosomes for degradation. Defects in this gene are a cause of susceptibility to inflammatory bowel disease type 10 (IBD10). Several transcript variants encoding different isoforms have been found for this gene.

抗原:Synthesized peptide derived from the Internal region of human Atg16

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: